

D-ACTOR® 100 »ultra«



D-ACTOR® mode:

- Radial acoustic wave:
pressure: 0.3 – 5.0 bar_{eff}
frequency: 1 – 21 Hz
(depending on pressure)
- Radial handpiece FALCON® with
integrated display and »Skin touch«
function for rapid treatment
navigation
- Vibration therapy (optional),
V-ACTOR®: 1 – 31 Hz
- 10" touch screen (optional)

D-ACTOR® 200 »ultra«



D-ACTOR® mode:

- Radial acoustic wave:
pressure: 0.3 – 5.0 bar_{eff}
frequency: 1 – 21 Hz
- Radial handpiece FALCON® with
integrated display and »Skin touch«
function for rapid treatment
navigation
- Two D-ACTOR® handpiece connectors
- Vibration therapy (optional),
V-ACTOR®: 1 – 35 Hz
- Vacuum therapy (optional),
VACU-ACTOR® technology
- 10" touch screen (optional)

CELLACTOR® SC1 T-TOP »ultra«



C-ACTOR® mode:

- Focused acoustic wave:
0.03 – 1.24 mJ/mm²
- Focused handpiece C-ACTOR®
SEPIA® with integrated display
for rapid treatment navigation
- 10" touch screen (optional)

CELLACTOR® SC1 TOWER »ultra«



C-ACTOR® mode:

- See CELLACTOR® SC1 T-TOP »ultra«

D-ACTOR® mode:

- See D-ACTOR® 200 »ultra«
- Integrated ultrasound imaging,
Colour Doppler (optional)
- Drawer module »Store Case«
(optional)
- 15,4" touch screen